

Sir Safety System S.p.A. Single-member company
Wholly-owned subsidiary of Sir Holding Srl
Via dei Fornaciai, 9 - Santa Maria degli Angeli
06081 Assisi - Perugia - Italy
Tel.+39 075.8043737 (R.A.) - Fax +39 075.8048686

VAT Nr. IT 03359340548
Viale Europa 74
20090 Cusago Milano
Tel. +39 02.90394575 (R.A.)
Fax +39 02.90394441

export@sirsafety.com
www.sirsafety.com

COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 9001 =



Name: Supplies for Medibag 2 With Sphygmomanometer	Code: 81149A	First Aid Equipment
---	---------------------	----------------------------

Description

Supplies kit for Medibag 2 with sphygmomanometer. First aid items required by D.M. 388 15/7/03 attachment 1 and DL 81 09/4/08 art.45.

Content:

- 1 copy D.M. 388 15/7/03
- 5 pairs sterile gloves
- 1 mask with face shield
- 3 sterile physiological solution, 500 ml
- 2 bottles antiseptic solution, 500 ml, 10% Povidone-iodine
- 10 packs of sterile compress 10x10 cm
- 2 packs sterile compress 18x40 cm
- 2 blankets 40x60 cm DIN 13152-BR
- 2 sterile forceps
- 1 tubular net gauze
- 1 pack cotton wool
- 2 Plastosan boxes of 10 assorted plasters
- 2 adhesive plaster rolls 5 m x 2.5 cm
- 1 pair of scissors 14.5 cm DIN 58279
- 3 haemostatic lacets
- 2 instant Ice Packs, single use
- 2 sanitary waste bags
- 1 digital clinical thermometer
- 1 Personal sphygmomanometer with sthethoscope
- 1 gauze bandage 3.5 m x 10 cm
- 1 First aid multilanguage instructions.



Category	Colours	Available Sizes Unica
Quantity/Carton 1	Important Information Read carefully the user sheet included. Make sure the item is in good condition, undamaged before use. Once you become aware of a damage, you should immediately replace it.	Recommended for Suitable for working environments, with three or more employees.
Dexterity	EC Standards	Storage and cleaning Keep in a dry place, far from heat sources.